
# European Gamma Sterilization and Aging Evaluations - Bromobutyls

# Small-Volume Serum Stoppers

- 13mm and 20mm
- Wide range of applications
- Elastomeric formulations meet requirements of major pharmacopoeia
- Elastomeric formulations for drugs with special characteristics
- Components can be siliconised and delivered ready-to-sterilise\*

\*If you intend to export your pharmaceutical products into the U.S. or Canada, West offers Westar® RS, a validated wash process with its own Type V DMF filed with the U.S. Food and Drug Administration and Health Canada.

Standard-grade rinsed product is available upon request.

Westar® is a registered trademark of West Pharmaceutical Services, Inc., in the United States and other jurisdictions.

# Small-Volume Lyophilization Stoppers

- 13mm and 20mm
- Wide range of applications
- Provide maximum stability for lyophilization to help ensure a safe and successful product
- Elastomeric formulations meet requirements of major pharmacopoeia
- Elastomeric formulations for drugs with special characteristics
- Components can be siliconised and delivered ready-to-sterilise\*

\*If you intend to export your pharmaceutical products into the U.S. or Canada, West offers Westar® RS, a validated wash process with its own Type V DMF filed with the U.S. Food and Drug Administration and Health Canada.

Standard-grade rinsed product is available upon request.

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## European Gamma Sterilization and Aging Evaluations - Bromobutyls

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<u>Formulation</u>	<u>Polymer</u>	<u>Cure System</u>
PH 4002/45 Gray	Bromobutyl	Resin
4023/50 Gray	Bromobutyl	Unconventional
4420/45 Light Gray	Bromobutyl	Unconventional
4405/50 Gray	Bromobutyl	Unconventional

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**European Gamma Sterilization and Aging Evaluations - Bromobutyls**

**Table I. Rubber Formulation PH 4002/45 Gray**

TESTS	SAMPLES			
	Control (No Processing)	Bierer Davies	25 kiloGray w/Bierer Davies	50 kiloGray w/Bierer Davies
<b>Tackiness</b> (Degree of Tackiness)	Medium	Strong	Strong	Strong
<b>Hardness</b> (Shore A Units)	48	46	44	45
<b>Microhardness</b> (IRHD)	52	45	43	42
<b>Compression Set</b> (22 hr/70°C/25% Deflection)	26	30	29	31
<b>Tensile Strength</b> (Mpa)	6.74	5.38	2.98	2.46
<b>Modulus</b> (100 Mpa)	1.10	0.78	0.78	0.91
<b>Elongation</b> (At Break %)	859	897	697	564
<b>EP Physical Properties</b>				
<b>Fragmentation</b> (Fragments/48 Punctures)	0	1	1	4
<b>Penetrability</b> (Newtons)	3.9	3.8	3.9	3.7
<b>EP Chemical Properties</b>				
<b>Appearance</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>pH Shift</b> (mL 0.01 NaOH/HCl)	0.04	0.04	0.04	0.04
<b>UV Absorption</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>Reducing Substances</b> (mL 0.01N KmnO <sub>4</sub> )	0.4	0.4	0.5	0.6
<b>Heavy Metals</b> (ppm)				
Pb	0	0	0	0
Zn	0.9	1.2	1.2	1.2
<b>Residue on Evaporation</b> (mg/50 mL)	0.0	0.0	0.0	0.0
<b>Volatile Sulphides</b> (mg Na <sub>2</sub> S/20 cm <sup>2</sup> )	0.00	0.00	0.00	0.00

**European Gamma Sterilization and Aging Evaluations - Bromobutyls**

<b>Table II. Rubber Formulation 4023/50 Gray</b>				
<b>TESTS</b>	<b>SAMPLES</b>			
	<b>Control (No Processing)</b>	<b>Bierer Davies</b>	<b>25 kiloGray w/Bierer Davies</b>	<b>50 kiloGray w/Bierer Davies</b>
<b>Tackiness</b> (Degree of Tackiness)	Slight	Slight	None	Slight
<b>Hardness</b> (Shore A Units)	50	50	55	54
<b>Microhardness</b> (IRHD)	51	51	50	51
<b>Compression Set</b> (22 hr/70°C/25% Deflection)	14	14	15	16
<b>Tensile Strength</b> (Mpa)	8.37	7.89	5.15	2.89
<b>Modulus</b> (100 Mpa)	2.77	3.10	2.91	3.11
<b>Elongation</b> (At Break %)	620	383	332	492
<b>EP Physical Properties</b>				
<b>Fragmentation</b> (Fragments/48 Punctures)	2	1	2	2
<b>Penetrability</b> (Newtons)	1.7	1.6	2.0	2.0
<b>EP Chemical Properties</b>				
<b>Appearance</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>pH Shift</b> (mL 0.01 NaOH/HCl)	0.0	0.0	0.0	0.0
<b>UV Absorption</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>Reducing Substances</b> (mL 0.01N KmnO <sub>4</sub> )	0.4	0.4	0.4	0.7
<b>Heavy Metals (ppm)</b>				
Pb	0	0	0	0
Zn	0.0	0.0	0.0	0.0
<b>Residue on Evaporation</b> (mg/50 mL)	0.0	0.0	0.0	0.0
<b>Volatile Sulphides</b> (mg Na <sub>2</sub> S/20 cm <sup>2</sup> )	0.06	0.06	0.00	0.00

**European Gamma Sterilization and Aging Evaluations - Bromobutyls**

**Table III. Rubber Formulation 4420/45 Light Gray**

TESTS	SAMPLES			
	Control (No Processing)	Bierer Davies	25 kiloGray w/Bierer Davies	50 kiloGray w/Bierer Davies
<b>Tackiness</b> (Degree of Tackiness)	None	Slight	Slight	Slight
<b>Hardness</b> (Shore A Units)	45	49	47	46
<b>Microhardness</b> (IRHD)	47	46	43	39
<b>Compression Set</b> (22 hr/70°C/25% Deflection)	18	31	25	38
<b>Tensile Strength</b> (Mpa)	7.19	6.92	5.80	3.88
<b>Modulus</b> (100 Mpa)	1.17	1.34	1.40	1.36
<b>Elongation</b> (At Break %)	787	737	755	639
<b>EP Physical Properties</b>				
<b>Fragmentation</b> (Fragments/48 Punctures)	0	1	1	1
<b>Penetrability</b> (Newtons)	1.5	2.0	1.7	1.5
<b>EP Chemical Properties</b>				
<b>Appearance</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>pH Shift</b> (mL 0.01 NaOH/HCl)	0.0	0.0	0.1	0.0
<b>UV Absorption</b>	Complies EP Type I	Complies EP Type I	Complies EP Type I	Complies EP Type I
<b>Reducing Substances</b> (mL 0.01N KmnO <sub>4</sub> )	0.1	0.4	0.0	0.2
<b>Heavy Metals</b> (ppm)				
Pb	0	0	0	0
Zn	0.0	0.0	0.0	0.0
<b>Residue on Evaporation</b> (mg/50 mL)	0.0	0.0	0.0	0.0
<b>Volatile Sulphides</b> (mg Na <sub>2</sub> S/20 cm <sup>2</sup> )	0.03	0.06	0.00	0.00

**European Gamma Sterilization and Aging Evaluations - Bromobutyls**

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- Standard washed product
- Enhanced performance product (EPP)
- Enhanced performance product, ready-to-sterilize (EPP RS)

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# Flip-Off® Seals

West's Flip-Off seals allow pharmaceutical manufacturers to run their filling machinery at high speed with less down time. Manufacturers using these high-precision seals may be able to reduce their pre-production inspection time.

## Features

- Tamper-evident closure
- Safe and easy to use with one-hand removal
- Suitable for high-speed processing
- Able to accommodate glass tolerances
- Helps to ensure closure integrity
- Steam-sterilisable
- Excellent cost/efficiency ratio
- Available in 13mm, 20mm, 28mm, 32mm

## Standard Flip-Off Button Colours

Black  
27607

Dark Gray  
4014

Natural (Clear)  
1000

White  
3766

Purple  
0527

Blue  
3769

Green  
3768

Yellow  
3772

Orange  
3773

Red  
3767

(Other colours are available upon request.)


# Quality of Injectables Starts with Quality of Primary Packaging Materials.

Exclusive Distributor in Thailand

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### BIOMED CO., LTD.