



VITREX® Soft BLOOD GLUCOSE LANCETS 23, 28, 30 and 32G

- ◆ RANGE OF BLOOD VOLUMES 23, 28, 30 and 32G lancets
- ◆ GOOD BLOOD FLOW triple ground needle point
- ◆ VERSATILE fit most commonly used lancing devices
- ◆ STERILE sterilised by gamma irradiation at minimum 25kGy
- ◆ CONVENIENT boxes have dispensing feature

PRODUCTS

DESCRIPTION	GAUGE SIZE	PIECES/BOX	ARTICLE No.
Sterile blood lancet	23G (0.65 mm)	200	351213
Sterile blood lancet	28G(0.38 mm)	200	352213
Sterile blood lancet	30G (0.31 mm)	200	354213
Sterile blood lancet	32G (0.26 mm)	200	357213

PACKAGING

PRIMARY BOX	CASE	SHIPPER	PALLET
200 pcs	10 x 200 2,000 pcs	10 x cases 20,000 pcs	480,000 pcs

PRODUCT SPECIFICATION

NEEDLE

MATERIAL	Stainless Steel Grade 304/1 Cr18Ni9		
EXPOSED LENGTH	3.3mm (+/- 0.33mm)		
POINT	3 facet (hypodermic) profile		
NOMINAL DIAMETER	0.65 mm (23G) 0.38 mm (28G) 0.31 mm (30G) 0.26 mm (32G)		

MOULDED BODY/CAP

MATERIAL	LDPE (LOW DENSITY POLYETHYLENE)		
OVERALL LENGTH	32.5 mm		
COLOUR	Light blue	23G	
	Red	28G	
	Green	30G	
	Dark blue	32G	

STERILITY

METHOD	Gamma irradiation
MINIMUM DOSAGE	25kGY
SHELF LIFE	5 years

VITREX® Soft IS DESIGNED FOR USE WITH MOST POPULAR TYPES OF FINGER LANCING DEVICE